

# **Bantam**

Petite but powerful 4F peripheral catheter 0.018" OTW featuring 6 fold on 7-9mm

## **NOW WITH 220MM LONG BALLOONS**

### INDICATIVE PRODUCT APPLICATIONS

The Bantam state-of-the-art 0.018" OTW PTA catheter designed for the optimum combination of profile, flexibility and pushability for the vast majority of upper and lower limb procedures.

#### **Durable Balloon**

The Bantam is equipped with ClearStream's acclaimed "QuadFlex" balloon noted for its durability in highly calcified vessels. The balloon is noted for its parallel surfaces, smooth inflation and perfect straightness and no "dog-bone".

#### **Longest Balloon**

The balloons are supplied in widest range of balloon sizes available on the market for PTA up to the world's longest at 220mm.

#### **Wide range of diameters**

The balloon has diameters ranging from 2.0mm to 9.0mm enabling the Bantam's unique features and performance to be brought to a broad range of procedures.

#### **Multiple shaft lengths for different approaches**

Bantam is supplied in a variety of shaft lengths to facilitate a wide range of access points and procedural approaches, including 75cm, 130cm and 180cm.

#### **Soft, low profile tip**

The Bantam's soft tip hugs the guidewire reducing vessel trauma and providing perfect trackability to deliver the Bantam's superior pushability.



The highlighted part numbers in the matrix indicate the Stocked Parts Inventory.

**Catheter Length 75cm**

Catheter Length	Balloon diameter	Balloon Length	Balloon Length	Balloon Length	Balloon Length	Balloon Length	Balloon Length	Balloon Length	Balloon Length
		20mm	40mm	60mm	80mm	100mm	120mm	150mm	220mm
75cm	2.0mm	<b>48402020</b>	48404020	48406020	48408020	48410020	48412020	48415020	48422020
75cm	2.5mm	48402025	48404025	<b>48406025</b>	48408025	48410025	48412025	48415025	<b>48422025</b>
75cm	3.0mm	48402030	48404030	<b>48406030</b>	48408030	48410030	48412030	<b>48415030</b>	48422030
75cm	3.5mm	48402035	<b>48404035</b>	48406035	<b>48408035</b>	<b>48410035</b>	48412035	48415035	48422035
75cm	4.0mm	48402040	48404040	48406040	48408040	<b>48410040</b>	<b>48412040</b>	48415040	<b>48422040</b>
75cm	5.0mm	48402050	<b>48404050</b>	48406050	<b>48408050</b>	<48410050	<b>48412050</b>	<b>48415050</b>	<b>48422050</b>
75cm	6.0mm	<b>48402060</b>	48404060	<b>48406060</b>	48408060	<b>48410060</b>	48412060	<b>48415060</b>	48422060
75cm	7.0mm		<b>48404070</b>	48406070					
75cm	8.0mm		<b>48404080</b>	48406080					
75cm	9.0mm		48404090	48406090					

### Catheter Length 90cm

Catheter Length	Balloon Diameter	Balloon Length	Balloon Length	Balloon Length	Balloon Length
		40mm	80mm	120mm	220mm
90cm	2.0mm	<b>48804020</b>	<b>48808020</b>	<b>48812020</b>	
90cm	2.5mm	<b>48804025</b>	<b>48808025</b>	<b>48812025</b>	<b>48822025</b>
90cm	3.0mm	<b>48804030</b>	<b>48808030</b>	<b>48812030</b>	<b>48822030</b>

### Catheter Length 130cm

Balloon Diameter	_____	_____	_____	Ballon	Length	_____	_____	_____	_____
	20mm	40mm	60mm	80mm	100mm	120mm	150mm	220mm	280mm
2.0mm	<b>48502020</b>	<b>48504020</b>	<b>48506020</b>	<b>48508020</b>	<b>48510020</b>	<b>48512020</b>	<b>48515020</b>	<b>48522020</b>	<b>48528020</b>
2.5mm	<b>48502025</b>	<b>48504025</b>	<b>48506025</b>	<b>48508025</b>	<b>48510025</b>	<b>48512025</b>	<b>48515025</b>	<b>48522025</b>	<b>48528025</b>
3.0mm	<b>48502030</b>	<b>48504030</b>	<b>48506030</b>	<b>48508030</b>	<b>48510030</b>	<b>48512030</b>	<b>4851030</b>	<b>48522030</b>	<b>48528030</b>
3.5mm	<b>48502035</b>	48504035	<b>48506035</b>	48508035	48510035	<b>48512035</b>	<b>48515035</b>	<b>48522035</b>	<b>48528035</b>
4.0mm	48502040	<b>48504040</b>	<b>48506040</b>	<b>48508040</b>	<b>48510040</b>	<b>48512040</b>	<b>48515040</b>	<b>48522040</b>	<b>48528040</b>
5.0mm	<b>48502050</b>	<b>48504050</b>	<b>48506050</b>	<b>48508050</b>	<b>48510050</b>	<b>48512050</b>	<b>48515050</b>	<b>48522050</b>	<b>48528050</b>
6.0mm	<b>48502060</b>	<b>48504060</b>	<b>48506060</b>	48508060	<b>48510060</b>	<b>48512060</b>	<b>48515060</b>	<b>48522060</b>	<b>48528060</b>
7.0mm		<b>48504070</b>	48506070						
8.0mm		48507080	48506080						
9.0mm		48504090	48506090						

### Catheter Length 150cm

Balloon Diameter	Balloon Length	Balloon Length	Balloon Length	Balloon Length
------------------	----------------	----------------	----------------	----------------

	40mm	80mm	100mm	120mm
5.0mm	47904050	47908050	47910050	47912050
6.0mm	47904060	47908060	47910060	47912060

**Catheter Length 180cm\***

Balloon Diameter	Balloon Length	Balloon Length	Balloon Length
	40mm	80mm	100mm
2.0mm	48004020	48008020	48010020
2.5mm	48004025	48008025	48010025
3.0mm	48004030	48008030	48010030
3.5mm	48004035	48008035	48010035
4.0mm	48004040	48008040	48010040

**\*180cm Length Catheter - special order only.**