

Two Dorado PTA Dilatation Catheters are shown diagonally against a light yellow background. The catheters are long, thin, and tapered towards the tip. The top one is slightly longer and more curved than the bottom one. Both have a dark, pointed tip.

# DORADO<sup>®</sup>

PTA Dilatation Catheter

Resistance Is Futile

# DORADO<sup>®</sup>

PTA Dilatation Catheter

## Confidence and precision for challenging cases

The DORADO<sup>®</sup> PTA Dilatation Catheter provides strength and versatility with enhanced deliverability, making this ultra non-compliant balloon the optimal choice for challenging lesions and in-stent restenosis.



CHECKER™ Flex Points

### ENHANCED DELIVERABILITY

- CHECKER™ Flex Points provide increased flexibility in tortuous anatomy
- Improved shaft design promotes faster inflation and deflation

### LONGER LENGTHS

For treatment of long stenosis or post-dilatation of long stents with fewer inflations.


Balloon Lengths to 20 cm

0 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20




## ULTRA NON-COMPLIANCE

- Concentrates maximum dilatation force at the resistant lesion
- Allows for inflation without the risk of overexpansion
- Strong balloon material promotes consistent post-stent dilatation
- Provides treatment of in-stent restenosis




Semi-compliant balloon in a resistant stenosis.



BARD® ultra non-compliant balloon in a resistant stenosis.

Balloon Compliance Comparison\*  
(7 mm Balloons)



\* Competitive compliance data taken from each manufacturer's product labeling.  
DORADO® PTA Dilatation Catheter compliance data on file.