


VALEO^{*}

VASCULAR STENT

powered to flexibility

STRENGTH



™ **BARD**

**PERIPHERAL
VASCULAR**



































Not for sale or distribution in the U.S.

VALEO*

VASCULAR STENT

Not for sale or distribution in the U.S.

Ordering Information

Expanded Stent Diameter [mm]	Mounted Stent Length [mm]	Product Code <i>for catheter working length of:</i>		Introducer Compatibility [F]	Nominal Pressure [atm]	RBP [atm]
		80 cm	120 cm			
6	18	EV06181C 	EV06182C 	6	8	14
	26	EV06261C 	EV06262C 	6	8	14
	36	EV06361C 	EV06362C 	6	8	14
7	18	EV07181C 	EV07182C 	6	8	14
	26	EV07261C 	EV07262C 	6	8	14
	36	EV07361C 	EV07362C 	6	8	14
	56	EV07561C 	EV07562C 	6	8	14
8	18	EV08181C 	EV08182C 	6	8	14
	26	EV08261C 	EV08262C 	6	8	14
	36	EV08361C 	EV08362C 	6	8	14
	56	EV08561C 	EV08562C 	7	8	14
9	17	EV09171C 	EV09172C 	7	6	12
	26	EV09261C 	EV09262C 	7	6	12
	36	EV09361C 	EV09362C 	7	6	12
10	17	EV10171C 	EV10172C 	7	6	12
	26	EV10261C 	EV10262C 	7	6	12
	36	EV10361C 	EV10362C 	7	6	12

For All Product Codes: 0.035" Guidewire.

PHYSICIAN'S SIGNATURE

AUSTRIA	Bard Medica SA, Thaliastrasse 125A/1/5, A-1160 Wien, Austria. Tel: +43 1 49 49130. Fax +43 1 49 49130 30.
BENELUX	Bard Benelux N.V., Hagelberg 2, B-2250 Olen, Belgium. Tel: +32 14 28 69 50. Fax +32 14 28 69 55.
FRANCE	Bard France S.A.S., Av Joseph Kessel 164-166, Parkile-P 14, BP 70, F-78960, Voisins-le-Bretonneux, France. Tel: +33 1 39 30 58 58. Fax +33 1 39 30 58 78.
GERMANY	C. R. Bard GmbH, Wachhausstrasse 6, D-76227 Karlsruhe, Germany. Tel: +49 721 9445 0. Fax: +49 721 9445 100.
GREECE	Bard Hellas SA, 22, Alkiviadou St & 72, Vouliagmenis Av, GR-16675 Glyfada, Greece. Tel: +30 210 9690770. Fax +30 210 9628810.
ITALY	Bard S.p.A, Via Cina 444, IT-00144 Roma, Italy. Tel: +39 06 52 493 222. Fax: + 39 06 52 95 454.
NORDIC	Bard Norden AB, Karbingatan 22, SE-254 67 Helsingborg, Sweden. Tel: +46 42 38 60 00. Fax: +46 42 38 60 10.
PORTUGAL	C. R. Bard Portugal Lda., R. Castilho, 13D-2 Piso A, 1250-066 Lisbon, Portugal. Tel: +351 213 190 330. Fax: +351 213 190 339.
SPAIN	Bard de Espana S.A., Avda Diagonal 652-656, Edificio A, 4 planta, E-08034 Barcelona, Spain. Tel: +34 93 253 78 00. Fax: +34 93 253 78 34.
SWITZERLAND	Bard Medica SA, Seestrasse 64, CH-8942 Oberrieden, Zurich, Switzerland. Tel: +41 44 722 53 60. Fax +41 44 722 53 70.
UK	Bard Limited, Forest House, Brighton Road, Crawley, West Sussex RH11 9BP, UK. Tel: +44 1293 52 95 55. Fax: +44 1293 52 83 98.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions, and information for use.

*Bard and Valeo are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.
© 2008 C. R. Bard, Inc. All Rights Reserved. S11649 Rev. 0



BARD

**PERIPHERAL
VASCULAR**